


STAMBP Antibody / AMSH / STAM binding protein [clone 28S66] (RQ8484)

Catalog No.	Formulation	Size
RQ8484	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28S66
Purity	Affinity chromatography
UniProt	O95630
Applications	Western Blot : 1:500-1:2000
Limitations	This STAMBP antibody is available for research use only.



Western blot testing of human HeLa cell lysate with STAMBP antibody. Predicted molecular weight ~48 kDa.

Description

STAM binding protein has a key role in regulation of cell surface receptor-mediated endocytosis and ubiquitin-dependent sorting of receptors to lysosomes. [UniProt]

Application Notes

Optimal dilution of the STAMBP antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human STAMBP protein was used as the immunogen for the STAMBP antibody.

Storage

Store the STAMBP antibody at -20°C.