


## Stap2 Antibody (R33749)

Catalog No.	Formulation	Size
R33749-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	106766
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Stap2 antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with Stap2 antibody at 1ug/ml. Predicted molecular weight: ~50kDa.

## Description

Additional name(s) for this target protein: Signal transducing adaptor family member 2

## Application Notes

Optimal dilution of the Stap2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids EEKLQKRRALEH were used as the immunogen for this Stap2 antibody.

## **Storage**

Aliquot and store the Stap2 antibody at -20oC.