


STAP2 Antibody (R33871)

Catalog No.	Formulation	Size
R33871-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55620
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This STAP2 antibody is available for research use only.



Western blot testing of human heart lysate with STAP2 antibody at 0.5ug/ml. Predicted molecular weight: ~50 kDa.

Description

Additional name(s) for this target protein: Signal transducing adaptor family member 2

Application Notes

Optimal dilution of the STAP2 antibody should be determined by the researcher.

1. This STP2 antibody is specific for isoforms 1 and 2 of the protein.

Immunogen

Amino acids ELQKKLEKRRALEH were used as the immunogen for this STAP2 antibody.

Storage

Aliquot and store the STAP2 antibody at -20oC.