


STAT5A Antibody (R33643)

Catalog No.	Formulation	Size
R33643-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6776
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This STAT5A antibody is available for research use only.



Western blot testing of human A) K562, B) MOLT4 and C) mouse NIH3T3 lysate with STAT5A antibody at 0.3ug/ml. Predicted molecular weight: ~91kDa.

Description

Additional name(s) for this target protein: Signal transducer and activator of transcription 5A

Application Notes

Optimal dilution of the STAT5A antibody should be determined by the researcher.

1. This STAT5A antibody will not cross react with STAT5B.

Immunogen

Amino acids DSRLSPPAGLFTSAR were used as the immunogen for this STAT5A antibody.

Storage

Aliquot and store the STAT5A antibody at -20oC.