


STAT5A Antibody (R34705)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34705-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow, Dog, Mouse, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 6776 |
| Applications | Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000 |
| Limitations | This STAT5A antibody is available for research use only. |



Western blot testing of K562 lysate with STAT5A antibody at 0.1ug/ml. Predicted molecular weight: ~91kDa.

Description

Additional name(s) for this target protein: Signal transducer and activator of transcription 5A

Application Notes

Optimal dilution of the STAT5A antibody should be determined by the researcher.

Immunogen

Amino acids KYYPVLAKAVD were used as the immunogen for this STAT5A antibody.

Storage

Aliquot and store the STAT5A antibody at -20oC.