


## STEAP4 Antibody (R33360)


Catalog No.	Formulation	Size
R33360-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79689
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This STEAP4 antibody is available for research use only.



Western blot testing of human adipose lysate with STEAP4 antibody at 0.3ug/ml.  
Predicted molecular weight: ~52kDa.



Western blot of HEK293 lysate overexpressing STEAP4 probed with STEAP4 antibody (mock transfection in lane 1). Predicted molecular weight: ~52kDa.

## Description

Additional name(s) for this target protein: STEAP family member 4; tumor necrosis factor, alpha-induced protein 9; TNFAIP9

## Application Notes

Optimal dilution of the STEAP4 antibody should be determined by the researcher.

## Immunogen

Amino acids CVDNTLTRIRQGWERN were used as the immunogen for this STEAP4 antibody.

## Storage

Aliquot and store the STEAP4 antibody at -20oC.