


## STK25 Antibody (F54066)

Catalog No.	Formulation	Size
F54066-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54066-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	O00506
<b>Applications</b>	Western Blot : 1:2000
<b>Limitations</b>	This STK25 antibody is available for research use only.



Western blot testing of STK25 antibody at 1:2000: Lane 1) human A431 cell lysate, 2) human HepG2 cell lysate and 3) human testis lysate. Predicted molecular weight ~48 kDa.

## Description

Oxidant stress-activated serine/threonine kinase that may play a role in the response to environmental stress. Targets to the Golgi apparatus where it appears to regulate protein transport events, cell adhesion, and polarity complexes important for cell migration. [UniProt]

## Application Notes

The stated application concentrations are suggested starting points. Titration of the STK25 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

Recombinant human STK25 was used as the immunogen for the STK25 antibody.

## **Storage**

Aliquot the STK25 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.