

STMN1 Antibody [Stathmin 1] (R34545)

Catalog No.	Formulation	Size
R34545-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3925
Applications	Western Blot : 1-3ug/ml Enzyme ImmunoAssay : binds recombinant protein ELISA (peptide) LOD : 1:1000
Limitations	This STMN1 antibody is available for research use only.



Western blot testing of human testis lysate with STMN1 antibody at 1ug/ml. Predicted molecular weight: ~17kDa.

Description

Additional name(s) for this target protein: Stathmin 1

Application Notes

Optimal dilution of the STMN1 antibody should be determined by the researcher.

Immunogen

Amino acids KNKESKDKPADETEAD were used as the immunogen for this STMN1 antibody.

Storage

Aliquot and store the STMN1 antibody at -20oC.