


STRA6 Antibody (R35189)


Catalog No.	Formulation	Size
R35189-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	20897
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This STRA6 antibody is available for research use only.



Western blot testing of siRNA-mediated STRA6 expression in 3T3-L1 lysate with STRA6 antibody at 1.7ug/ml.



Western blot: The 90kDa band is stained after retrovirally overexpressing mouse STRA6 and not after overexpressing human STRA6 in 3T3 preadipocytes.

Description

Additional name(s) for this target protein: Stimulated by retinoic acid gene 6

Application Notes

Optimal dilution of the STRA6 antibody should be determined by the researcher.

Immunogen

Amino acids QAFRKAALTSAKAN were used as the immunogen for this STRA6 antibody.

Storage

Aliquot and store the STRA6 antibody at -20oC.