


STX1A/B Antibody (R34383)

Catalog No.	Formulation	Size
R34383-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6804; 112755
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This STX1A/B antibody is available for research use only.



Western blot testing of mouse brain lysate with STX1A/B antibody at 0.003ug/ml.
Predicted molecular weight: ~33kDa (STX1A and STX1B).

Description

Additional name(s) for this target protein: Syntaxin 1A; Syntaxin 1B; STX1B

Application Notes

Optimal dilution of the STX1A/B antibody should be determined by the researcher.

Immunogen

Amino acids KDRTQELRTAKD were used as the immunogen for this STX1A/B antibody. This sequence is common to STX1A and STX1B.

Storage

Aliquot and store the STX1A/B antibody at -20oC.