


STXBP6 Antibody (R35832)

Catalog No.	Formulation	Size
R35832-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29091
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This STXBP6 antibody is available for research use only.



IHC testing of FFPE human kidney with STXBP6 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Syntaxin binding protein 6, Amisyn

Application Notes

Optimal dilution of the STXBP6 antibody should be determined by the researcher.

Immunogen

Amino acids SAKSAISKEIFAP were used as the immunogen for this STXBP6 antibody.

Storage

Aliquot and store the STXBP6 antibody at -20oC.