


SWAP70 Antibody (R33760)

Catalog No.	Formulation	Size
R33760-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23075
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SWAP70 antibody is available for research use only.



Western blot testing of human lymph node lysate with SWAP70 antibody at 0.1ug/ml.
Predicted molecular weight: ~70kDa.

Description

Additional name(s) for this target protein: Switching B-cell complex 70kDa

Application Notes

Optimal dilution of the SWAP70 antibody should be determined by the researcher.

Immunogen

Amino acids EEREKNWKEKKTTE were used as the immunogen for this SWAP70 antibody.

Storage

Aliquot and store the SWAP70 antibody at -20oC.