


## Synaptogyrin 4 Antibody (R34205)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34205-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                                  |
|--------------------|------------------------------------------------------------------|
| Availability       | 1-3 business days                                                |
| Species Reactivity | Human                                                            |
| Format             | Antigen affinity purified                                        |
| Host               | Goat                                                             |
| Clonality          | Polyclonal (goat origin)                                         |
| Isotype            | Goat Ig                                                          |
| Purity             | Antigen affinity                                                 |
| Gene ID            | 23546                                                            |
| Applications       | Western Blot : 0.01-0.1ug/ml<br>ELISA (peptide) LOD : 1:128000   |
| Limitations        | This Synaptogyrin 4 antibody is available for research use only. |



Western blot testing of human brain lysate with Synaptogyrin 4 antibody at 0.01ug/ml.  
Predicted molecular weight: ~26kDa.

## Description

Additional name(s) for this target protein: SYNGR4

## Application Notes

Optimal dilution of the Synaptogyrin 4 antibody should be determined by the researcher.

## Immunogen

Amino acids HIPKSLQELANS were used as the immunogen for this Synaptogyrin 4 antibody.

## **Storage**

Aliquot and store the Synaptogyrin 4 antibody at -20oC.