


Synaptotagmin-2 Antibody (R34562)

Catalog No.	Formulation	Size
R34562-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	127833
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Synaptotagmin-2 antibody is available for research use only.



Western blot testing of mouse lung lysate with Synaptotagmin-2 antibody at 0.1ug/ml.
Predicted molecular weight: ~47kDa.

Description

Additional name(s) for this target protein: SYT2

Application Notes

Optimal dilution of the Synaptotagmin-2 antibody should be determined by the researcher.

Immunogen

Amino acids KPEEEVDALLGKNK were used as the immunogen for this Synaptotagmin-2 antibody.

Storage

Aliquot and store the Synaptotagmin-2 antibody at -20oC.