


# Synaptotagmin-2 Antibody (R34562)

Catalog No.	Formulation	Size
R34562-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	127833
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Synaptotagmin-2 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: SYT2

## **Application Notes**

Optimal dilution of the Synaptotagmin-2 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids KPEEEVDALLGKNK were used as the immunogen for this Synaptotagmin-2 antibody. **Storage** Aliquot and store the Synaptotagmin-2 antibody at -20oC.