


## Synaptotagmin-5 Antibody (R36112)

Catalog No.	Formulation	Size
R36112-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6861
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Synaptotagmin-5 antibody is available for research use only.



Western blot testing of human brain (olfactory bulb) lysate with Synaptotagmin-5 antibody at 1ug/ml. Predicted molecular weight: ~43kDa.

## Description

Additional name(s) for this target protein: Syt5, synaptotagmin V, SytV; Syt5

## Application Notes

Optimal dilution of the Synaptotagmin-5 antibody should be determined by the researcher.

## Immunogen

Amino acids SYIDKVQPEIEELD were used as the immunogen for this Synaptotagmin-5 antibody.

## **Storage**

Aliquot and store the Synaptotagmin-5 antibody at -20oC.