


SYNGAP1 Antibody (R33984)

Catalog No.	Formulation	Size
R33984-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8831
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SYNGAP1 antibody is available for research use only.



Western blot testing of mouse brain lysate with SYNGAP1 antibody at 0.03ug/ml.
Predicted/observed molecular weight: ~148kDa.

Description

Additional name(s) for this target protein: Synaptic Ras GTPase activating protein 1

Application Notes

Optimal dilution of the SYNGAP1 antibody should be determined by the researcher.

Immunogen

Amino acids ESAHIEREEYKLKEY were used as the immunogen for this SYNGAP1 antibody.

Storage

Aliquot and store the SYNGAP1 antibody at -20oC.