


Syntaxin 6 Antibody (R35196)


Catalog No.	Formulation	Size
R35196-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10228
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Syntaxin 6 antibody is available for research use only.



Western blot testing of human lymph node lysate with Syntaxin 6 antibody at 0.3ug/ml.
Predicted molecular weight: ~29kDa.



Western blot of HEK293 lysate overexpressing human STX6 probed with Syntaxin 6 antibody (mock transfection in lane 1). Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: STX6

Application Notes

Optimal dilution of the Syntaxin 6 antibody should be determined by the researcher.

Immunogen

Amino acids QALAERKNRQA were used as the immunogen for this Syntaxin 6 antibody.

Storage

Aliquot and store the Syntaxin 6 antibody at -20oC.