


Syntenin Antibody (R35956)

Catalog No.	Formulation	Size
R35956-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6386
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Syntenin antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with Syntenin antibody at 0.01ug/ml.
Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: SDCBP, Syndecan binding protein

Application Notes

Optimal dilution of the Syntenin antibody should be determined by the researcher.

1. This Syntenin antibody will detect isoforms 1 and 3.

Immunogen

Amino acids SLEDLKVDKVIQAQ were used as the immunogen for this Syntenin antibody.

Storage

Aliquot and store the Syntenin antibody at -20oC.