


Syntenin Antibody (R35956)

Catalog No.	Formulation	Size
R35956-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6386
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Syntenin antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with Syntenin antibody at 0.01ug/ml.
Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: SDCBP, Syndecan binding protein

Application Notes

Optimal dilution of the Syntenin antibody should be determined by the researcher.

1. This Syntenin antibody will detect isoforms 1 and 3.

Immunogen

Amino acids SLEDLKVDKVIQAQ were used as the immunogen for this Syntenin antibody.

Storage

Aliquot and store the Syntenin antibody at -20oC.