


TACR1 Antibody (R33806)

Catalog No.	Formulation	Size
R33806-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6869
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TACR1 antibody is available for research use only.



Western blot testing of human testis lysate with TACR1 antibody at 1ug/ml. Predicted molecular weight ~46kDa but routinely observed at ~58kDa. Image shows detection of the 'short' isoform at ~36kDa.

Description

Additional name(s) for this target protein: Tachykinin receptor 1

Application Notes

Optimal dilution of the TACR1 antibody should be determined by the researcher.

Immunogen

Amino acids EIPGDSSDRYHEQ were used as the immunogen for this TACR1 antibody.

Storage

Aliquot and store the TACR1 antibody at -20oC.