


TACR1 Antibody (R33806)

Catalog No.	Formulation	Size
R33806-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6869
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TACR1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Tachykinin receptor 1

Application Notes

Optimal dilution of the TACR1 antibody should be determined by the researcher.

Immunogen

Amino acids EIPGDSSDRYHEQ were used as the immunogen for this TACR1 antibody. **Storage** Aliquot and store the TACR1 antibody at -20oC.