


## TAIP-12 Antibody (R34632)

Catalog No.	Formulation	Size
R34632-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81566
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TAIP-12 antibody is available for research use only.



Western blot testing of human heart lysate with TAIP-12 antibody at 0.1ug/ml. Predicted molecular weight: ~59kDa.

## Description

Additional name(s) for this target protein: TGF-beta induced apoptosis protein 12, FAM130A1; BC035295

## Application Notes

Optimal dilution of the TAIP-12 antibody should be determined by the researcher.

## Immunogen

Amino acids KRQKQLRRKNVRFDQ were used as the immunogen for this TAIP-12 antibody.

## Storage

Aliquot and store the TAIP-12 antibody at -20oC.