


## Tankyrase 2 Antibody (R34273)

Catalog No.	Formulation	Size
R34273-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80351
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Tankyrase 2 antibody is available for research use only.



Western blot testing of Jurkat lysate with Tankyrase 2 antibody at 0.5ug/ml. The expected ~125kDa band and the additional ~45kDa band are both blocked by the immunizing peptide.

### Description

Additional name(s) for this target protein: TRF1-interacting ankyrin-related ADP-ribose polymerase 2; TNKS2

### Application Notes

Optimal dilution of the Tankyrase 2 antibody should be determined by the researcher.

## Immunogen

Amino acids HRRKEVSEENHNH were used as the immunogen for this Tankyrase 2 antibody.

## Storage

Aliquot and store the Tankyrase 2 antibody at -20oC.