


## Tankyrase 2 Antibody (R34273)

Catalog No.	Formulation	Size
R34273-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80351
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Tankyrase 2 antibody is available for research use only.



Western blot testing of Jurkat lysate with Tankyrase 2 antibody at 0.5ug/ml. The expected ~125kDa band and the additional ~45kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: TRF1-interacting ankyrin-related ADP-ribose polymerase 2; TNKS2

## Application Notes

Optimal dilution of the Tankyrase 2 antibody should be determined by the researcher.

## Immunogen

Amino acids HRRKEVSEENHHH were used as the immunogen for this Tankyrase 2 antibody.

## **Storage**

Aliquot and store the Tankyrase 2 antibody at -20oC.