


Tankyrase 2 Antibody (R35581)

Catalog No.	Formulation	Size
R35581-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80351
Applications	Western Blot : 1-3ug/ml
Limitations	This Tankyrase 2 antibody is available for research use only.



Western blot testing of transfected HEK293 lysate with Tankyrase 2 antibody at 0.5ug/ml. Lane 1: mock TANK2; Lane 2: TANK1; Lane 3: empty; Lane 4: TANK2. Lower panel shows the same lysates probed for Actin (loading control).

Description

Additional name(s) for this target protein: TRF1-interacting ankyrin-related ADP-ribose polymerase 2, TNKS2

Application Notes

Optimal dilution of the Tankyrase 2 antibody should be determined by the researcher.

1. This Tankyrase 2 antibody will not detect Tankyrase 1.

Immunogen

A portion of amino acids 800-900 was used as the immunogen for this Tankyrase 2 antibody.

Storage

Aliquot and store the Tankyrase 2 antibody at -20oC.