


Tankyrase Antibody (R35580)

Catalog No.	Formulation	Size
R35580-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8658
Applications	Western Blot : 0.5-1.5ug/ml
Limitations	This Tankyrase antibody is available for research use only.



Western blot of HEK293 lysate overexpressing TANK1/Tankyrase (lane 2) and TANK2 (lane 4) tested with Tankyrase antibody at 0.5ug/ml. Mock transfection is in lane 1 and lane 3 is empty. Lower panel shows the same lysates tested for alpha-Actin (loading control).

Description

Additional name(s) for this target protein: TRF1-interacting ankyrin-related ADP-ribose polymerase, TNKS

Application Notes

Optimal dilution of the Tankyrase antibody should be determined by the researcher.

Immunogen

A recombinant protein fragment of 120 amino acids from the N-terminus was used as the immunogen for this Tankyrase

antibody.

Storage

Aliquot and store the Tankyrase antibody at -20oC.