


TAZ Antibody (R34229)

Catalog No.	Formulation	Size
R34229-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6901
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TAZ antibody is available for research use only.



Description

Additional name(s) for this target protein: Tafazzin

Application Notes

Optimal dilution of the TAZ antibody should be determined by the researcher.

1. This TAZ antibody will detect isoforms 1/2/3/4 of the protein.

Immunogen

Amino acids HLKTQAEQLHNH were used as the immunogen for this TAZ antibody.

Storage

Aliquot and store the TAZ antibody at -20oC.