

## TAZ Antibody (R34229)

Catalog No.	Formulation	Size
R34229-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6901
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TAZ antibody is available for research use only.



Western blot testing of human heart lysate with TAZ antibody at 0.3ug/ml. Predicted molecular weight: 28~33kDa (isoforms 1-4).

## Description

Additional name(s) for this target protein: Tafazzin

## Application Notes

Optimal dilution of the TAZ antibody should be determined by the researcher.

1. This TAZ antibody will detect isoforms 1/2/3/4 of the protein.

## **Immunogen**

Amino acids HLKTQAEQLHNH were used as the immunogen for this TAZ antibody.

## **Storage**

Aliquot and store the TAZ antibody at -20oC.