


TBC1D3E Antibody (F40022)

Catalog No.	Formulation	Size
F40022-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40022-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	A6NER0
Applications	Western Blot : 1:1000
Limitations	This TBC1D3E antibody is available for research use only.



Description

TBC1D3E acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11 (By similarity).

Application Notes

Titration of the TBC1D3E antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 498-527 from the human protein was used as the immunogen for this TBC1D3E antibody.

Storage

Aliquot the TBC1D3E antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.