


TBK1 Antibody (R36186)

Catalog No.	Formulation	Size
R36186-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29110
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TBK1 antibody is available for research use only.



Western blot testing of HeLa lysate with TBK1 antibody at 1ug/ml. The expected ~85kDa band and the additional ~24kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: TANK-binding kinase 1; NAK; NF-kB-activating kinase

Application Notes

Optimal dilution of the TBK1 antibody should be determined by the researcher.

Immunogen

Amino acids TIETSLQDIDSRLS were used as the immunogen for this TBK1 antibody.

Storage

Aliquot and store the TBK1 antibody at -20oC.