


TBLR1 Antibody (R35842)

Catalog No.	Formulation	Size
R35842-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6907; 79718
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This TBLR1 antibody is available for research use only.



IHC testing of FFPE human kidney with TBLR1 antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Transducin (beta)-like 1X-linked receptor 1, TBL1X, TBL1XR1

Application Notes

Optimal dilution of the TBLR1 antibody should be determined by the researcher.

Immunogen

Amino acids SASDGSVCVLDLRK were used as the immunogen for this TBLR1 antibody.

Storage

Aliquot and store the TBLR1 antibody at -20oC.