


TBP Antibody (R33598)

Catalog No.	Formulation	Size
R33598-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6908
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TBP antibody is available for research use only.



Western blot testing of HeLa nuclear [A] and cytosolic [B] lysate with TBP antibody at 0.3ug/ml. Predicted molecular weight ~38kDa.

Description

Additional name(s) for this target protein: TATA box binding protein

Application Notes

Optimal dilution of the TBP antibody should be determined by the researcher.

1. This TBP antibody is specific for isoform 1 of the protein.

Immunogen

Amino acids DQNN~~S~~LPPYAQ were used as the immunogen for this TBP antibody.

Storage

Aliquot and store the TBP antibody at -20°C.