


TCF19 Antibody (R35780)

Catalog No.	Formulation	Size
R35780-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6941
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TCF19 antibody is available for research use only.



Western blot testing of 293 lysate with TCF19 antibody at 1.5ug/ml. Predicted molecular weight: ~39kDa.

Description

Additional name(s) for this target protein: Transcription factor 19, SC1

Application Notes

Optimal dilution of the TCF19 antibody should be determined by the researcher.

Immunogen

Amino acids RSTAKAPSDTPAHE were used as the immunogen for this TCF19 antibody (accession # AAB60363.1).

Storage

Aliquot and store the TCF19 antibody at -20oC.