


## TCF3 Antibody / Transcription factor E2-alpha Biotin Conjugate (R34355)

Catalog No.	Formulation	Size
R34355BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Biotin Conjugate
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	6929
<b>Applications</b>	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This TCF3 antibody is available for research use only.



Western blot testing of Daudi cell lysate with biotinylated TCF3 antibody at 0.3ug/ml.  
Expected molecular weight: 68-78 kDa.

### Description

Additional name(s) for this target protein: Transcription factor 3, ITF1

### Application Notes

Optimal dilution of the TCF3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids KAPRARTSPDEDED were used as the immunogen for this TCF3 antibody.

## **Storage**

Aliquot and store the TCF3 antibody at -20oC.