


TCF4 Antibody (R34798)

Catalog No.	Formulation	Size
R34798-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6925
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TCF4 antibody is available for research use only.



Western blot testing of human heart lysate with TCF4 antibody at 0.5ug/ml. Predicted molecular weight ~71kDa.

Description

Additional name(s) for this target protein: Transcription factor 4, E2-2, ITF2

Application Notes

Optimal dilution of the TCF4 antibody should be determined by the researcher.

1. This TCF4 antibody will detect isoforms a and b.

Immunogen

Amino acids NNGHPSPSRNYGDGT were used as the immunogen for this TCF4 antibody.

Storage

Aliquot and store the TCF4 antibody at -20°C.