


TCFL5 Antibody (R34701)

Catalog No.	Formulation	Size
R34701-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10732
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TCFL5 antibody is available for research use only.



Western blot testing of human testis lysate with TCFL5 antibody at 0.3ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: Transcription factor-like 5

Application Notes

Optimal dilution of the TCFL5 antibody should be determined by the researcher.

Immunogen

Amino acids KYIQERHGDSLKKE were used as the immunogen for this TCFL5 antibody.

Storage

Aliquot and store the TCFL5 antibody at -20oC.