


TCHHL1 Antibody (R35197)

Catalog No.	Formulation	Size
R35197-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	126637
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This TCHHL1 antibody is available for research use only.



Western blot testing of human skin lysate with TCHHL1 antibody at 1ug/ml. Predicted molecular weight: ~99kDa.

Description

Additional name(s) for this target protein: Trichohyalin-like 1

Application Notes

Optimal dilution of the TCHHL1 antibody should be determined by the researcher.

Immunogen

Amino acids QALEDKQGHPQRE were used as the immunogen for this TCHHL1 antibody.

Storage

Aliquot and store the TCHHL1 antibody at -20oC.