


TCL1 Antibody (R35369)

Catalog No.	Formulation	Size
R35369-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8115
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TCL1 antibody is available for research use only.



Western blot testing of Jurkat cell lysate with TCL1 antibody at 0.1ug/ml. Predicted molecular weight: ~14kDa.

Description

Additional name(s) for this target protein: T-cell leukemia/lymphoma 1A; TCL1A

Application Notes

Optimal dilution of the TCL1 antibody should be determined by the researcher.

Immunogen

Amino acids QLYPDGRYRSSD were used as the immunogen for this TCL1 antibody.

Storage

Aliquot and store the TCL1 antibody at -20oC.