


# TCL1 Antibody (R35369)

Catalog No.	Formulation	Size
R35369-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8115
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TCL1 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: T-cell leukemia/lymphoma 1A; TCL1A

### **Application Notes**

Optimal dilution of the TCL1 antibody should be determined by the researcher.

### **Immunogen**

Amino acids QLYPDGRYRSSD were used as the immunogen for this TCL1 antibody.

# **Storage** Aliquot and store the TCL1 antibody at -20oC.