


TERF1 Antibody (R34854)

Catalog No.	Formulation	Size
R34854-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7013
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TERF1 antibody is available for research use only.



Western blot testing of human testis lysate with TERF1 antibody at 0.3ug/ml. Predicted molecular weight: ~50/65kDa (unmodified/possible acetylated form).

Description

Additional name(s) for this target protein: Telomeric repeat binding factor 1, TRF1, PIN2

Application Notes

Optimal dilution of the TERF1 antibody should be determined by the researcher.

1. This TERF1 antibody is specific for isoform 1.

Immunogen

Amino acids NLDTRKSVSDKQS were used as the immunogen for this TERF1 antibody.

Storage

Aliquot and store the TERF1 antibody at -20°C.