


TFB1M Antibody (R33866)

Catalog No.	Formulation	Size
R33866-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51106
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TFB1M antibody is available for research use only.



Western blot testing of A431 lysate with TFB1M antibody at 0.3ug/ml. Predicted molecular weight: ~39kDa.

Description

Additional name(s) for this target protein: Transcription factor B1, mitochondrial

Application Notes

Optimal dilution of the TFB1M antibody should be determined by the researcher.

Immunogen

Amino acids ELKRRKSKNEEKE were used as the immunogen for this TFB1M antibody.

Storage

Aliquot and store the TFB1M antibody at -20oC.