


TFB2M Antibody (R34457)

Catalog No.	Formulation	Size
R34457-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64216
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TFB2M antibody is available for research use only.



Western blot testing of human lung lysate with TFB2M antibody at 0.1ug/ml.
Predicted/observed molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Transcription factor B2, mitochondrial

Application Notes

Optimal dilution of the TFB2M antibody should be determined by the researcher.

Immunogen

Amino acids KHCFGRRSATVIDH were used as the immunogen for this TFB2M antibody.

Storage

Aliquot and store the TFB2M antibody at -20°C.