


## TFEB Antibody (R35316)

Catalog No.	Formulation	Size
R35316-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7942
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This TFEB antibody is available for research use only.



Western blot testing of human lung lysate with TFEB antibody at 1ug/ml. Expected molecular weight: 53-60 kDa (unmodified), 65-70 kDa (modified).

## Description

Additional name(s) for this target protein: Transcription factor EB

## Application Notes

Optimal dilution of the TFEB antibody should be determined by the researcher.

## **Immunogen**

Amino acids QKDLQKSRELENH were used as the immunogen for this TFEB antibody.

## **Storage**

Aliquot and store the TFEB antibody at -20oC.