


TFF3 Antibody / Trefoil factor 3 (R35836)

Catalog No.	Formulation	Size
R35836-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7033
Localization	Cytoplasmic, secreted
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TFF3 antibody is available for research use only.



Western blot testing of colon cancer lysate with TFF3 antibody at 1ug/ml. Predicted molecular weight: ~10 kDa.

Description

Additional name(s) for this target protein: Trefoil factor 3

Application Notes

Optimal dilution of the TFF3 antibody should be determined by the researcher.

Immunogen

Amino acids SANQCAVPAKDRVD were used as the immunogen for this TFF3 antibody.

Storage

Aliquot and store the TFF3 antibody at -20oC.