


TFPI Antibody / Tissue factor pathway inhibitor (R36331)


Catalog No.	Formulation	Size
R36331-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7035
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TFPI antibody is available for research use only.



Western blot testing of HeLa lysate with TFPI antibody at 0.3ug/ml. Predicted/observed molecular weight: 35/40-50kDa.



Western blot of HEK293 lysate overexpressing TFPI and probed with TFPI antibody (mock transfection in lane 1). Predicted/observed molecular weight: 35/40-50kDa.

Description

Additional name(s) for this target protein: Tissue Factor Pathway Inhibitor

Application Notes

Optimal dilution of the TFPI antibody should be determined by the researcher.

1. This TFPI antibody is specific for isoform a.

Immunogen

Amino acids VKIAYEEIFVKNM were used as the immunogen for this TFPI antibody.

Storage

Aliquot and store the TFPI antibody at -20oC.