


TGFBI Antibody (R35372)

Catalog No.	Formulation	Size
R35372-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7045
Applications	Western Blot : 0.2-0.5ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TGFBI antibody is available for research use only.

kDa 200 —	Western blot testing of human kidney lysate with TGFBI antibody at 0.3ug/ml. Predicted
116 – 97 –	molecular weight: ~75 kDa.
66 -	
44	
31 –	
22 –	
14 –	
6-	



IHC staining of FFPE human kidney with TGFBI antibody at 3.8 μ ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Transforming growth factor, beta-induced, 68kDa

Application Notes

Optimal dilution of the TGFBI antibody should be determined by the researcher.

Immunogen

Amino acids QLYTDRTEKLRPE were used as the immunogen for this TGFBI antibody.

Storage

Aliquot and store the TGFBI antibody at -20oC.