


TGIF2 Antibody (R33882)

Catalog No.	Formulation	Size
R33882-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	60436
Applications	Western Blot : 0.3-1.0ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TGIF2 antibody is available for research use only.



Western blot testing of A431 lysate with TGIF2 antibody at 0.3ug/ml. Predicted molecular weight: ~26kDa.

Description

Additional name(s) for this target protein: TGFB-induced factor 2

Application Notes

Optimal dilution of the TGIF2 antibody should be determined by the researcher.

Immunogen

Amino acids EMELQKQDP were used as the immunogen for this TGIF2 antibody.

Storage

Aliquot and store the TGIF2 antibody at -20oC.