


THBS3 Antibody / Thrombospondin-3 (R33284)

Catalog No.	Formulation	Size
R33284-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7059
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This THBS3 antibody is available for research use only.



IHC testing of FFPE human colon with THBS3 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Thrombospondin 3; TSP3

Application Notes

Optimal dilution of the THBS3 antibody should be determined by the researcher.

Immunogen

Amino acids CNDTVPEDFEPFRRQ were used as the immunogen for this THBS3 antibody.

Storage

Aliquot and store the THBS3 antibody at -20oC.