


THNSL2 Antibody (R33781)

Catalog No.	Formulation	Size
R33781-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55258
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This THNSL2 antibody is available for research use only.



Western blot testing of human colorectal cancer lysate with THNSL2 antibody at 2ug/ml.
Predicted molecular weight: ~45kDa.

Description

Additional name(s) for this target protein: Threonine synthase-like 2

Application Notes

Optimal dilution of the THNSL2 antibody should be determined by the researcher.

Immunogen

Amino acids EGNDELDEPIKT were used as the immunogen for this THNSL2 antibody.

Storage

Aliquot and store the THNSL2 antibody at -20oC.